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## User Manual

## INDEX

<b>1</b>	<b>GENERAL</b> .....	<b>3</b>
1.1	COPYRIGHT .....	3
1.2	WARNING AND CAUTIONING SIGNS .....	3
1.3	DISCLAIMER .....	3
1.4	MANUFACTURER .....	4
1.5	CE MARK, APPLICABLE RULES AND REGULATIONS.....	4
1.6	USER INSTRUCTIONS FOR SAFE OPERATION AND OPERATOR QUALIFICATION REQUIREMENTS .....	4
<b>2</b>	<b>USEFUL LIFE AND WARRANTY CONDITIONS</b> .....	<b>5</b>
<b>3</b>	<b>SCOPE OF DELIVERY</b> .....	<b>5</b>
<b>4</b>	<b>INTENDED USE</b> .....	<b>5</b>
<b>5</b>	<b>MOUNTING AND COMMISSIONING</b> .....	<b>6</b>
<b>6</b>	<b>ELECTRICAL CONNECTIONS</b> .....	<b>6</b>
<b>7</b>	<b>MOUNTING</b> .....	<b>6</b>
<b>8</b>	<b>PRODUCT DESCRIPTION AND CONTROLS</b> .....	<b>7</b>
<b>9</b>	<b>OPERATION OF CHAIR</b> .....	<b>8</b>
1.1.	BRAKES .....	8
1.2.	BACKREST .....	8
1.2.1.	ARMRESTS .....	8
1.2.2.	HEADREST .....	8
1.3.	SEAT .....	8
1.4.	HEIGHT ADJUSTENT (OPTIONAL).....	8
<b>10</b>	<b>CLEANING AND CONTAMINATION PROTECTION</b> .....	<b>9</b>
<b>11</b>	<b>MAINTENANCE</b> .....	<b>9</b>
<b>12</b>	<b>SAFETY INSPECTIONS</b> .....	<b>9</b>
<b>13</b>	<b>DISPOSAL OF DEVICE</b> .....	<b>9</b>
<b>14</b>	<b>TECHNICAL DATA</b> .....	<b>10</b>
<b>15</b>	<b>ENVIRONMENTAL CONDITIONS</b> .....	<b>10</b>

## 1 General

To safely operate and control the medical device please carefully read this user manual before operating. Please pay attention to all cautionary statements pertaining to the safe operation of the medical device. Please keep the manual on file at all times while the product is in service.

The manufacturer reserves the right to make technical changes or modifications in accordance with advanced technical developments. This manual is not subject to an updating service.

### 1.1 Copyright

© Distributing, copying or any commercial use of this manual is strictly prohibited unless expressly permitted in writing by the manufacturer. The manufacturer is entitled to compensation for any violation of this right.

The manufacturer reserves any and all rights in case a patent is granted or the industrial design is registered.

### 1.2 Warning and cautioning signs

	Warning	Indicates a hazardous situation which may result in fatal or serious bodily injury if the appropriate safety precautions are not heeded.
	Caution	Indicates a hazardous situation which may result in minor injury if the appropriate safety precautions are not heeded.  Caution Property damage Indicates possible device damage if the appropriate safety precautions are not heeded
		This pictograph on the label indicates Device typ B according to EN 60601-1

### 1.3 Disclaimer

	Warning	This device may only be set up, operated and used for the specified purpose and according to national regulations, consistent with the applicable industry standards and occupational safety and accident prevention regulations..
---	---------	--

This declaration shall be rendered invalid if changes are made to the product without the manufacturer's authorization.

1.4 Manufacturer

**AKRUS GmbH & Co KG**

Otto-Hahn-Straße 3

D-25337 ELSHORN

☎ int. +49 (0) 4121 791930

Email: [info@akrus.de](mailto:info@akrus.de)

FAX int. +49 (0) 4121 791939

Website: [www.akrus.de](http://www.akrus.de)

1.5 CE mark, applicable rules and regulations

The device complies with the following Regulations for medical devices.

- MDD 93/42 EEC
- DIN / ISO 60601-1
- DIN / ISO 60601-1-2 EMC
- UMDNS No 13-961
- Medical Device Class I

The device meets the requirements of MDD 93/42/EEC and the corresponding national requirements of the Medical Devices Act MDA (Medizinprodukte Gesetz MPG).

User must observe the national legal requirements governing the prevention of accidents.

1.6 User instructions for safe operation and operator qualification requirements

**Before using this medical device, please carefully read and observe the safety instructions and recommendations listed in this manual**

	<p>Caution Risk of operating errors</p>	<p>The device may only be installed, operated, used and maintained by persons who have been properly trained or who have the required knowledge and experience to do so. Please also adhere to the national qualification guidelines applicable in your country.</p>
---	---	--

- Proper handling and use of the medical device is essential to its safe operation.
- Operation of the medical device only in accordance with the intended use.
- Do not operate the medical device in explosive or hazardous environments. The medical device must not be operated, with the exception of very minute quantities, in the presence of combustible anesthetics, or volatile solvents like alcohol, benzene or the like.
- Do not set up the medical device in humid or damp rooms. Never expose the product to sprinkling, dripping or splashing water.
- Pay special attention to and carefully read any emphasized safety advice or information in this user manual.

## 2 Useful life and warranty conditions

The useful product life is expected to be 8 years. Both, the useful life time and the warranty offered by manufacturer are conditional upon the adherence to the conditions outlined below.

This declaration shall be rendered invalid if changes are made to the product without the manufacturer's authorization.

	<p>Caution Risk of operating errors</p>	<p>The device may only be installed, operated, used and maintained by persons who have been properly trained or who have the required knowledge and experience to do so. Please also adhere to the national qualification guidelines applicable in your country.</p>
	<p>Warning</p>	<p>The development, production and maintenance of these instrument tables, together with associated risks, are based on an expected service life of eight years, provided the device is serviced at the specified intervals. Modifications to the product or failure to follow the manufacturer's instructions may substantially reduce the expected service life and significantly increase the risks associated with the use of this device. It is the responsibility of the institution operating this product to follow the manufacturer's instructions and to decide on the risk/benefit ratio upon expiration of the expected service life or maintenance and inspection intervals specified by the manufacturer.</p>
	<p>Warning General Hazards</p>	<p>This device may only be set up, operated and used for the specified purpose and according to national regulations, consistent with the applicable industry standards and occupational safety and accident prevention regulations. Further notes on classification are to be found in chapter <i>Technical data</i>.</p>

## 3 Scope of delivery

The ak 5003M scope of delivery comprises

Position	Qty
ak 5003M mammography chair	1
Power cord	1
Quality- test report	1
User Manual	1
Options according to order	

## 4 Intended use

The medical device is intended for use as a surgical chair. Any other use of the product is not permitted.

- The medical device is designed for a maximum patient weight of 150 kg (330 lbs) (EN 60601-2-46 particular requirements for safety of operating tables). The headrest is designed for a maximum load of 20 kg (44 lbs).
- The chair is suitable to be employed in the immediate proximity of patients.

## 5 Mounting and commissioning

The medical device is delivered completely assembled. If optional parts or accessories are part of the delivery, please refer to the respective section in this user manual for mounting or connecting instructions.

If not expressly mentioned in this user manual, no additional assembly or calibration work is required to operate the medical device.

## 6 Electrical connections

The location of the power input of the device or the recharging unit is described in this user manual.

To connect the medical device to the power supply use only the power cable included in this delivery.

	<p>Warning General Hazard</p>	<p>The medical device is only de-energized when the mains plug has been pulled.</p> <p>The device should be set up so that the power cable can be disconnected from the power supply quickly and easily without any tools.</p>
	<p>Warning Risk of electric shock</p>	<p>Additional extension cables or portable multiple sockets may not be used.</p> <p>The electrical installation must conform to IEC 60364-7-710 or the applicable national regulations. This includes the existence of a ground fault circuit interrupter (GFCI).</p> <p>To avert the risk of an electric shock this device may only be connected to a power supply with a protective earth conductor.</p> <p>Ensure that the power supply plug is suitable and certified for the local connection. If the supplied power cable must be replaced, at least the following specifications must be adhered to:</p> <ul style="list-style-type: none"> <li>• Protective ground conductor resistance of the power cable maximum 0.1 Ohm</li> <li>• Local certification of the power cable for connection to medical devices</li> <li>• Device plug C19 conforming to IEC 60320</li> <li>• Cross-section at least 1.5 mm<sup>2</sup>/AWG 16</li> </ul>
	<p>Warning Fire Hazard</p>	<p>The device is not suitable for operation in explosion risk areas (e.g. combustible mixture of anesthetic, cleaning or disinfecting agents with air, oxygen or nitrous oxide).</p> <p>The electrical installation must conform to IEC 60364-7-710.</p> <p>The data concerning power consumption on the type plate must be taken into account in the selection of over current protection.</p>

## 7 Mounting

After mounting connect powercord with the chair. The connector is under the baseplate. Put plug into socket

	<p>Caution</p>	<p>Don't roll over the powercord. Look after the cord. Use only cords in good condition. Don't move chair without disconnecting.</p>
---	----------------	--

## 8 Product description and controls

The patient chair ak 5003M is to be used for examinations in the Mammography- area.

The ak 5003M has a castor under every corner of the baseplate. Every castor has its own brake.


There is a knob on both sides of the backrest for adjustment. The backrest is manually reclinable, down to its horizontal position. The footrest moves with the backrest. It converts to enlarge the lying area.


The headrest is infinitely adjustable. When not used it is detachable.

The armrests move along with the reclination and they are moveable.

Under the seat is a lever on both sides to turn around the seat in 360°. The seat is infinitely adjustable.

Optional: There are 2 footswitches on the baseplate for height adjustment.

## 9 Operation of chair


Treat a patient on this chair only when all 4 castors are fixed.

### 1.1. Brakes

Push down the brakelever for fixing the castor. To release push the lever upwards.

### 1.2. Backrest

Recline the backrest by pushing down the knob. Release the knob and the backrest is fixed.

#### 1.2.1. Armrests


For easy entry move the armrests up.

#### 1.2.2. Headrest

Move the headrest up or down for height adjustment. To detach it pull the cord out of the backrest. For attachment push the weight at the end of the cord into the backrest gently.

### 1.3. Seat

To turn around the seat loosen the lever under the seat. After turning tighten the lever for fixation.

### 1.4. Height adjustent (Optional)

Control the height adjustment of the seat via footswitch; the left switch stands for up and the right one for going down

## 10 Cleaning and contamination protection

- Contamination protection  
In order to protect the product against contamination by liquids, body liquids or other unwanted substances it is recommended to, while in use, cover the product with a non sterile disposable impermeable cover sheet.
- Surfaces are resistant to all commonly used disinfectants as listed in the D G H M (Deutsche Gesellschaft fuer Hygiene und Mikrobiologie) list.  
<http://www.gesundheitsamt.de/alle/technik/download/seuche/bakt/list/dghm.pdf>. Any appropriate disinfectant must not exceed the concentrations given for the following ingredients. Propanol=35% // Ethanol=25%.
- Do not attempt to sterilize the medical device.
- Clean only the external surfaces of the equipment using a damp cloth. In case of heavy stains do **not** use abrasive or aggressive materials other than commonly used cleansers or detergents. Do not use ether, acetone or concentrated acids.
- No water or cleaning detergents may enter the inside of the equipment.
- To clean soiled upholstery surfaces, use a damp (not wet) cloth and/or commonly used upholstery spray.
- Before resuming operation, allow sufficient time to dry.

## 11 Maintenance

The medical device is maintenance-free to users and operators.

The product is subject to the regular service intervals. Refer to the current issue of manufacturers "maintenance record".

## 12 Safety inspections

No safety inspections are mandated by the manufactures for this medical device. However, refer to the local (national) "Minimum requirements for safety inspections for medical devices class I".

## 13 Disposal of device

In accordance with applicable EU guidelines and national regulations at the time at which the product was brought onto the market, the product specified on the consignment note is not to be disposed of via the domestic waste disposal system or communal waste disposal facilities.

For further information on disposing of this product please contact your local dealer or the manufacturer or its legal successor company. Please read the latest information provided on the manufacturer's website.

If the product or its components are resold, the seller must inform the buyer that the product must be disposed of in accordance with the currently applicable national regulations.

#### 14 Technical data

Technical Data	Value	Unit
<b>Dimensions and weight</b>		
Length wheel base (backrest up)	700	mm
Width wheel base	645	mm
Width seat cushion	590	mm
Height chair back rest upright	1515	mm
Max. length chair, back rest horizontal	1865	mm
Max. patient weight	150	kg
Chair weight (depending on options)	96	kg
<b>Vertical range chair (Z axis)</b>		
Lowest Entry height	605	mm
Lift	200	mm
<b>Range back rest</b>		
Vertical to horizontal	0-90	°
<b>Wheel base</b>		
4 wheels		
- wheels pivoting separately		
- wheels locked separately		
Wheel diameter	100	mm
<b>Electrical data</b>		
Line power 230V/120V	230/ 120	V
Nominal frequency 230V/120V	50/ 60	Hz
Nominal power	3	A
Fuse	T5/20 3, 15A	A/H//V
Protection category	IP 65	
Continuous operation of electric motors	1	min
Upholstery electric conductible (optional)		

#### 15 Environmental Conditions

<b>Ambient conditions for intended use</b>	
Ambient air temperature	+10° C to +40° C
Rel. Humidity max	50 % non-condensing
Ambient air pressure	700 – 1060 hPa
<b>Ambient conditions for storage</b>	
Ambient air temperature	-10 °C to +55 °C
Rel. Humidity	10 % to 95 % non-condensing
<b>Ambient conditions for storage and transport in original packaging</b>	
Ambient air temperature	-40 °C to +70 °C
Rel. Humidity	10 % to 95 % non-condensing

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