Arzthocker Doctor Stool

ak 445 / ak 446 / ak 447

ak 447

ak 446

ak 445

Gebrauchsanweisung User Manual

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User Manual

Index

1. GENERAL.....	3
1.1. COPYRIGHT.....	3
1.2. WARNING AND CAUTIONING SIGNS.....	3
1.3. DISCLAIMER.....	3
1.4. MANUFACTURER	4
1.5. CE DECLARATION, APPLICABLE RULES AND REGULATIONS	4
1.6. USER INSTRUCTIONS FOR SAFE OPERATION	4
2. USEFUL LIFE AND WARRANTY CONDITIONS	5
3. SCOPE OF DELIVERY.....	5
4. INTENDED USE	6
5. MOUNTING AND COMMISSIONING	6
6. PRODUCT DESCRIPTION AND CONTROLS.....	6
6.1. HEIGHT ADJUSTMENT.....	6
6.2. BACK REST	6
6.3. ARM RESTS	6
7. OPERATION OF DOCTOR STOOL.....	7
7.1. BASE AND CASTORS.....	7
7.2. HEIGHT ADJUSTMENT.....	7
7.2.1. HEIGHT ADJUSTMENT, HAND OPERATED	7
7.2.2. HEIGHT ADJUSTMENT, FOOT OPERATED (AK 445 & AK 447).....	7
7.2.3. HEIGHT ADJUSTMENT, FOOT OPERATED (AK 446)	8
7.3. ADJUSTMENT OF SEAT TILT (AK 447 ONLY).....	8
7.4. BACK REST	8
7.4.1. ADJUSTMENT OF TILT ANGLE	8
7.4.2. ADJUSTMENT OF HEIGHT	9
7.5. ARM RESTS	9
7.5.1. STANDARD ARM REST	9
7.5.2. MULTI ARTICULATED ARM REST	10
7.5.2.1. LENGTH EXTENSION OF ARM RESTS.....	10
7.5.2.2. TILT OF ARM REST	10
7.5.2.3. HEIGHT ADJUSTMENT OF ARM REST	10
8. CLEANING AND ROUTINE CARE	11
9. MAINTENANCE.....	11
10. SAFETY INSPECTIONS	11
11. DISPOSAL OF DEVICE	11
12. TECHNICAL DATA.....	12

1. General

To safely operate and control this medical product please carefully read this user manual before operating the chair. Please pay attention to all cautionary statements pertaining to the safe operation of this medical product. Please keep the manual on file at all times while the product is in service.

The manufacturer reserves the right to make technical changes or modifications in accordance with advanced technical developments. This manual is not subject to an updating service.

1.1. Copyright

© Distributing, copying or any commercial use of this manual is strictly prohibited unless expressly permitted in writing by the manufacturer. The manufacturer is entitled to compensation for any violation of this right.

The manufacturer claims any and all rights in case a patent is granted or the design is registered.

1.2. Warning and cautioning signs

	Warning	Indicates a hazardous situation which may result in fatal or serious bodily injury if the appropriate safety precautions are not heeded.
	Caution	<p>Indicates a hazardous situation which may result in minor injury if the appropriate safety precautions are not heeded.</p> <p>Caution Property damage Indicates possible device damage if the appropriate safety precautions are not heeded</p>

1.3. Disclaimer

	Warning	This device may only be set up, operated and used for the specified purpose and according to national regulations, consistent with the applicable industry standards and occupational safety and accident prevention regulations..
---	---------	--

Manufacturer excludes any and all liabilities in case of improper or not authorized operation or maintenance or servicing of the product.

This declaration shall be rendered invalid if changes are made to the product without the manufacturer's authorization.

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☎ int. +49 (0) 4121 791930
 FAX int. +49 (0) 4121 791939

Email: info@akrus.de
 Website: www.akrus.de

1.5. CE declaration, applicable rules and regulations

The device complies with the following Regulations for machine products.

- 2006/42/EG
- RL 2011/65/ EU

The device meets the requirements of 2006/42/EG and the corresponding national requirements of the machinery directive.

User must observe the national legal requirements governing the prevention of accidents.

1.6. User instructions for safe operation

Before using this medical product, please carefully read and observe the safety instructions and recommendations listed in this manual

	Caution Risk of operating errors	The device may only be installed, operated, used and maintained by persons who have been properly trained or who have the required knowledge and experience to do so. Please also adhere to the national qualification guidelines applicable in your country.
---	---	---

- To safely use and operate the medical device, use only as directed.
- Operation of the medical device only in accordance with the intended use.
- Do not operate the medical device in explosive or hazardous environments. The medical device must not be operated, with the exception of very minute quantities, in the presence of combustible anesthetics, or volatile solvents like alcohol, benzene or the like.
- Do not set up the medical device in humid or damp rooms. Never expose the product to sprinkling, dripping or splashing water.
- Modifications, authorized by manufacturer, may only be performed by staff authorized by manufacturer. Any modification non authorized or instructed by manufacturer may cause malfunctions or hazards to persons.
- Maintenance or repair work on this product may only be performed by staff authorized by manufacturer.
- Pay special attention to any emphasized safety instruction or information in this user manual.
- Only trained and authorized staff may operate the medical device.
- Only trained service staff authorized by manufacturer may open the medical device.

2. Useful life and warranty conditions

The useful product life is expected to be 8 years. Both, the useful life time and the warranty offered by manufacturer are conditional upon the adherence to the conditions outlined below.

This declaration shall be rendered invalid if changes are made to the product without the manufacturer's authorization.

	<p>Caution Risk of operating errors</p>	<p>The device may only be installed, operated, used and maintained by persons who have been properly trained or who have the required knowledge and experience to do so. Please also adhere to the national qualification guidelines applicable in your country.</p>
---	---	--

	<p>Warning</p>	<p>The development, production and maintenance of these instrument tables, together with associated risks, are based on an expected service life of eight years, provided the device is serviced at the specified intervals. Modifications to the product or failure to follow the manufacturer's instructions may substantially reduce the expected service life and significantly increase the risks associated with the use of this device. It is the responsibility of the institution operating this product to follow the manufacturer's instructions and to decide on the risk/benefit ratio upon expiration of the expected service life or maintenance and inspection intervals specified by the manufacturer.</p>
---	----------------	---

	<p>Warning General Hazards</p>	<p>This device may only be set up, operated and used for the specified purpose and according to national regulations, consistent with the applicable industry standards and occupational safety and accident prevention regulations. Further notes on classification are to be found in chapter <i>Technical data</i>.</p>
---	--	--

3. Scope of delivery

The ak 445/ 446/ 447 scope of delivery comprises

Position	No
Doctor stool ak 445, ak 446 or ak 447	1
Quality-test-report	1
User manual	1
Options according to order	

4. Intended use

The ak 445/ 446/ 447 doctor stools are designed to offer maximum ergonomic comfort for medical procedures performed in a seated position. The stools are not approved for the transport of people. Any other use of the product is not permitted.

- The stools are designed for a maximum load of 135 kg.
- The chair is suitable to be used in the immediate proximity of patients.

5. Mounting and commissioning

The medical device is delivered completely assembled. If optional parts or accessories are part of the delivery, please refer to the respective section in this user manual for mounting or connecting instructions.

If not expressly mentioned in this user manual, no additional assembly or calibration work is required to operate the medical device.

6. Product description and controls

The ak 445 and 447 doctor stools are available with either a round seat cushion (ak 445, ak 446) or a saddle seat (ak 447). All adjustments are manually.

6.1. Height adjustment

The doctor stools are available with either hand or foot operated height adjustment.

6.2. Back rest

The optional back rest of the doctor stools features one lever for height adjustment and one lever to lock the spring loaded in a fixed position or release it to follow the operators movements.

6.3. Arm rests

Arm rests are available for the ak 445 and ak 446 models. Available are either standard or multi articulated arm rests.

7. Operation of doctor stool

7.1. Base and Castors

The doctor stools are equipped with 5 electrically conductible heavy duty castors allowing effortless movements of the stool. A “heel kick” brake is available as an option.

	<p>Warning Gen. hazard</p>	<p>To move the chair over a door step or threshold ≥ 20 mm, PULL chair at the back rest, thus avoiding any danger of tipping.</p>
---	------------------------------------	---

7.2. Height adjustment

The stool is height adjustable by a pneumatic or hydraulic spring. The spring is unlocked by either a hand operated lever or foot operated (heel kick) stainless steel plate (ak 445 & ak 447).

After unlocking the gas spring (hand or foot operated) the body weight will push down the stool. To move the stool up, remove the body weight off the stool (ak 445 & ak 447).

The stool moves up by pushing the foot lever. Push the lever down to end, the hydraulic spring will be released and the stool moves downwards (ak 446).

The stool will continue to move up or down until either the control is set to „Lock“ or the stool is in full stop position (ak 445 & ak 447).

7.2.1. Height adjustment, hand operated

The hand operated lever for both the ak 445 and 447 is located under the seat cushion. Pull up and hold until desired seating position was reached.

7.2.2. Height adjustment, foot operated (ak 445 & ak 447)

The foot operated release plate is located underneath the base. Apply gentle pressure by the heel to unlock the gas spring. Press and hold down the plate until desired seating position was reached.

	<p>Warning Gen. hazard</p>	<p>Never kick down with both heels on two different spots. This way the gas spring cannot be unlocked and also, this may cause damage to the mechanism.</p>
---	------------------------------------	---

7.2.3.Height adjustment, foot operated (ak 446)

Kick down the heel kick lever to operate the hydraulic pump to elevate the chair up. To lower the chair, press and hold down the lever beyond the soft stop barrier. (ak 446).

Heel kick lever

7.3. Adjustment of seat tilt (ak 447 only)

The seat cushion of the ak 447 doctor stool offers a tilt of approximately 10°.The lever is located underneath the seat.

Tilt max. bwd

Lever

Tilt max. fwd

7.4.1.Adjustment of tilt angle

The release lever for the back rest tilt is located under the seat. If pulled into the upper (horizontal) position, the spring loaded back rest will follow the movements of the back rest. To lock the back rest into any desired position, pushed lever to the lower position (see picture).

Lever

7.4.2. Adjustment of height

The height of the back rest is adjustable using the lever at the base of the back rest.

- Horizontal position, back rest slides up - down
- Vertical position, back rest is locked

7.5. Arm rests

Different arm rests are available for the ak 445 and ak 446, which can be adjusted to the ergonomic demands of the operator. The lower supporting arms swivel forward and back ward.

Swiveling brackets

7.5.1. Standard arm rest

The standard arm rests are height adjustable and, in addition to the swivel of the lower supporting arm, pivot in and out.

7.5.2. Multi articulated arm rest

The multi articulated arm rests offer unlimited possibilities to adjust the arm rest to any ergonomically desired or otherwise necessary position, thus preventing fatigue and offering relaxing working conditions. The adjustments are controlled by easily accessible locking handles.

	Warning	Make sure to properly lock arm rest in position to avoid any unwanted movements.
---	---------	--

7.5.2.1. Length extension of arm rests

To lock or unlock the upper part of arm rest use the locking handle at the side. Arm rest slides forward or back ward by approximately 20 cm.

7.5.2.2. Tilt of arm rest

To tilt arm rest up – down or sideways, use the locking handle at the base of the ball head.

7.5.2.3. Height adjustment of arm rest

For height adjustments of the arm rests or swivel of the upper part use the locking handle at the back of the lower supporting arm.

8. Cleaning and routine care

- Contamination protection

In order to protect the product against contamination by liquids, body liquids or other unwanted substances it is recommended to, while in use, cover the product with a non sterile disposable impermeable cover sheet.

- Surfaces are resistant to all commonly used disinfectants as listed in the D G H M (Deutsche Gesellschaft für Hygiene und Mikrobiologie) list. www.dghm.org
- Do not attempt to sterilize the medical device.
- Clean only the external surfaces of the equipment using a damp cloth. In case of heavy stains do **not** use abrasive or aggressive materials other than commonly used cleansers or detergents. Do not use ether, acetone or concentrated acids.
- No water or cleaning detergents may enter the inside of the equipment.
- To clean soiled upholstery surfaces, use a damp (not wet) cloth and/or commonly used upholstery spray.
- Before resuming operation, allow sufficient time to dry.

9. Maintenance

The chair is maintenance-free to users and operators.

Product is subject to the regular service intervals. Refer to the current issue of local (national) “Minimum requirements for medical products class I”.

10. Safety Inspections

No safety inspections are mandated by the manufactures for this medical device. However, refer to the local (national) “Minimum requirements for safety inspections for medical devices class I”.

11. Disposal of device

All components are rated house hold waste.

12. Technical data

Technical Data	ak 445	ak 447	ak 446	
Diameter base	600	600	600	mm
Diameter castors	50	50	50	mm
Height adjustment	530 - 720	530 - 720	550 - 750	mm
Seat cushion	420	430 x 340	420	mm
Thickness seat cushion	80 mm	80 mm	80	mm
Height adjustment arm rest	140 - 290	140 - 290	140 - 290	mm
Max. body weight	135	135	135	kg
Environmental conditions In Use				
Ambient air temperature	+10 - +40	+10 - +40	+10 - +40	° C
Relative air humidity max	50	50	50	%
Ambient air pressure	700 - 1060	700 - 1060	700 - 1060	hPa