LSC-Config

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Inhalt

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Figure Index

<i>Fig. 1: LSC-Config. Summary Tab</i>	2
<i>Figure 2: Settings</i>	3
<i>Figure 3: Status bar</i>	3
<i>Figure 4: Tab Joystick</i>	4
<i>Figure 5: Adjustment of joystick</i>	4
<i>Figure 6: Tab Axis Data</i>	5
<i>Figure 7: Tab Diagnosis</i>	6
<i>Figure 8: Tab Options</i>	6
<i>Figure 9: Example of a Summary-file</i>	7

1 Intended Use

The software is used to capture and visualize data related to the design, production and/or service for the LSC patient bed to be used with a refractive laser.

2 General Description


Fig. 1: LSC-Config. Summary Tab

The application covers essentially 4 areas (see. fig.1):

1. Menu bar (top)
2. Tabs bar (left)
3. Data-/control area (center)
4. Status bar (bottom)

The menu bar comprises File-, Settings-, Connection- and About-menu with additional menu sub items. The tabs bar includes the following tabs: summary, joystick, axis data, diagnostics and options. The controls for bed settings, input- and display elements for data are contained in the data/ control area. Information about the status of the connection to and the state of the bed are being displayed in the status bar.

The connection can be established either through the icon *Connect*  , in the lower part of the tab bar, or through the menu > Connection.

2.1 Setting up A CONNECTION

To set up a connection to the patient bed key in the settings first. Go to the menu bar *Settings* -> *Setting*..


Figure 2: Settings

In the window *Settings* (see figure. 2) choose the communication port (*COM interface*), corresponding transfer rate (*Baud rate*) and a time window to receive a message. Confirm the settings by clicking the icon *Set*. These settings are now stored and automatically uploaded when using the device. D

Figure 3: Status bar

Now click the icon *Connect* . After the connection is established the indicator in the status bar shows „*Connection status:*“ denomination of the COM-Ports (see Fig.. 3) and on the far right of the status bar it shows *Active*. The icon *Connect* changes to *Disconnect* .

3 Tabs

The controls for bed settings and the input and display elements for all data are located in the tabs.

3.1 Tab Joystick

The tab (see Fig. 4) contains the data and settings for the joystick as well as the status indication for each icon and the inverse-LED. The zero, min and max. position for each axis as well as axis deflection as a percentage indication are shown here. Additionally, in this tab a zero-spot window (zero range) can be set for each axis. A joystick adjustment is also available in this tab.


Figure 4: Tab Joystick

Adjustment of joysticks

To start the joystick adjustment, press the *Start*-Icon in the drop-down menu *Options*. After pressing the icon, a popup menu appears showing the instructions for the joystick adjustment. There are two phases of adjustment: zero position and max./min position.


Figure 5: Adjustment of joystick

The adjustment of the joystick can be stopped at any phase by pressing CANCEL.

3.2 Tab Axis Data


Figure 6: Tab Axis Data

The tab *Axis Data* (fig.. 6) contains all data and settings for all axes. For the X-, Y- and Z-axis these positions are displayed numerical and graphically. For each axis the status indications for the current status of this axis (*Enabled*) and the validity of data (*Valid*) as well as the current speed, PWM-value and the power consumption are being displayed in this tab. The max. travel distance limitations and the PWM values can be altered in their respective input boxes.

The current status, the lock position and the power consumption of the locking actuator are indicated in this tab.

3.3 Tab Diagnosis

The individual POST-steps (POST Overview) as well as the entire POST-progress are being displayed in the tab *Diagnosis*. The system status indications like POST Supervisor- and Interlock-Status are also located here.

Any error in POST it can be tracked POST Overview.


Figure 7: Tab Diagnosis

3.4 Tab Options


Figure 8: Tab Options

The tab Options comprises:

- Automatic Functions
- Interlock Settings
- Default Settings
- Service Settings

The automatic functions can be set in *Automatic Functions*. Choose one or more axes for each execution phase and set an individual target position for each axis.

A read only access is available for the *Interlock Settings* in the service mode. The max. PWM values and travel distances with Interlock 1 in the active mode are being displayed here.

To return to factory settings go to icon *Factory Defaults* in *Default Settings*.

Enter the serial number and date of manufacturing in the *Service Settings*.

4 Save Data

To save the general data of the bed and the modules from the data files tab *Summary* (fig.. 1) as a PDF file, go to File -> Save. To enter the path where to save the file go to Settings -> Settings... -> *Storage path*.

Data of LSC80

Serial Number: -----
 Manufacturing Date: --:--

Modules

Modul	Serial Number	Software Version
M380	-----	745
UM	-----	688
UM	-----	688
UM	-----	688
APM2	-----	643
APM2	-----	643
APM3	-----	643

Figure 9: Example of a Summary-file

5 Known Errors

- Indication *Not connected* in the Status bar: a connection cannot be established, the r COM-Port is possibly busy.
- Missing data in the tab *Summary*: the bed was possibly switched on after it being connected to LSC Config -> re connect Software