



Product information
ak445 / ak446 / ak447



Doctor's stools



ak445




ak446

ak447

Pneumatic height adjustment by foot or hand release

Saddle seat, tiltable

AKRUS stool upholstery is of manufactory quality



Here you will find the technical data for the doctor's stools

Doctor's stools

The classic doctor's stool, the pump stool and the saddle stool

AKRUS doctor's stools are versatile and flexible medical seating furnitures. They are intended to support healthcare professionals during patient examination or patient treatment.

The classic among doctor's stools is the **ak445**. Equipped with a round seat, the height is adjusted using a gas pressure spring, which can be ordered either with foot or hand release. Optional, comfortable and flexible armrests can be easily adjusted to the desired position. Horizontally movable armrest holders allow you to navigate the armrests with your elbows.

If a pump stool is needed, the pump stool **ak446** comes into play. Its height is adjusted using a hydraulic pump system. A 360 degree pumping lever allows for ease of use even during a medical procedure.

The saddle stool **ak447** provides an ergonomically pre-shaped saddle like seat cushion. As with the doctor's stool ak445, the saddle stool ak447 is height adjusted using a gas pressure spring, which can be ordered either with foot or hand release. It can be completed with backrest and armrests.

All stools are moved by a stable roller cross holding five high-quality rollers. They can optionally be equipped with pre-shaped, multi-directional, multi-joint armrests. This allows even finer tuning of the forearm/hand support during an operation. The AKRUS upholstery is of manufactory quality. Various colors are available to match the color of the clinic or practice environment.

The doctor's stools ak445, ak446 and ak447 are delivered fully assembled and the use of high-quality components allows for many years of reliable use when used as intended. In the event of a service requirement, authorized medical technology dealers are available.



Multi articulated armrests



Foot release



Saddle seat cushion

Focus: ergonomic comfortable sitting

ak445



Classic doctor's stool



Foot release

ak446



Pump stool, front view



Pump stool, back view

ak447



Saddle Stool, front view



Saddle Stool, back view

ak445



Backrest adjustment levers



Plain stool with hand release

ak446



Multi articulated armrest with various fine tuning options



360° foot pump lever

ak447



Hand release levers for height and tilt adjustments

Upholstery fabrics

Choose one of 16 colors from our standard range. Special colors on request.



Combra White /
207x605




Creme /
208x2380




Sand /
207x3261



Pearl /
032x4354



Lilac /
208x2392



Bluegrey /
208x1277



Corn flower /
207x729



Midnightblue /
032x4709



Yellow /
208x2795



Orange /
208x3107




Red /
207x2796




Oper /
218x1637




Mint /
208x2394



Turquoise /
208x749



Charcoal /
207x1049



Black /
electrically conductive
anti-static
217x3117

Distributor